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Tel.: 03731 392835

[Juliane.Grahl@zuv.tu-freiberg.de](mailto:Juliane.Grahl@zuv.tu-freiberg.de)

Institute of Thermal-, Environmental- and Resources' Process Engineering

## Enhancement of Raman signals in aqueous samples

### BACKGROUND

Raman spectroscopy enables contactless chemical analysis of aqueous samples, for example in process monitoring in the chemical industry or in food and medical analytics. However, because Raman scattering is inherently weak, the resulting signals are often difficult to interpret. Signal-enhancing technologies are therefore essential for improving analytical performance.

### TECHNICAL DESCRIPTION


Fig. 2: SEM image of a silica aerogel coating inside an ALC

The core of the technology is a measurement cell consisting of a glass capillary whose inner wall is coated with a water-repellent silica aerogel (Fig. 1). The highly porous coating of such an aerogel-lined capillary (ALC) has an air-like refractive index (i.e. lower than that of aqueous samples), enabling light guiding through the liquid sample inside the capillary via total internal reflection. Both the excitation light and the Raman scattering generated within the sample are guided back to the Raman sensor.

In contrast to classical cuvette measurements, which evaluate only a tiny focal spot, ALCs capture the entire sample volume — resulting in a substantial signal enhancement. In addition, the hydrophobic aerogel coating protects the capillary from corrosion and contamination.

### APPLICATION AREA

Raman spectroscopy of aqueous samples

### ADVANTAGES

- ✓ enhancement of Raman signals in aqueous samples
- ✓ only small sample volumes required (minimum approx. 0.5 µL)
- ✓ no interference signals from the measurement cell itself

### STATUS

- ✓ proof of functionality demonstrated using industrial aqueous samples
- ✓ invention disclosure submitted

### COLLABORATION OPTIONS

- ✓ patent licensing
- ✓ R&D cooperation
- ✓ contract research

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